


GOAL: Client will be able to turn toothbrush on/off while it is in his mouth (to prevent spraying saliva and toothpaste while bristles are rotating).

Client has Braun Oral B toothbrush similar to this one. Here are the dimensions of his toothbrush.



The client has significant extensor tone throughout his extremities. When using precision or cylindrical grasps, the client tends to flex the distal joints of his fingers but keeps the proximal joints straight. This results in less contact overall with the object. See photos on next page.

With time and effort, he can turn the toothbrush on/off when he is holding it in front of him. He cannot turn it on/off when it is in his mouth.

